

BD® LSR II Flow Cytometer, Core Unit Fluoreszenz-Technologien, Universität Leipzig

Laser (nm)	Filter	Fluorochrome	Alternatives
Violet (405)	450/50	BV421	Pacific Blue
	525/50	BV510	AmCyan, VioGreen
Blue (488)	488/10	SSC	
	530/30	FITC	AF488, CFSE, GFP
	575/26	PE	PI
	695/40	PerCP-Cy5.5	
	780/60	PE-Cy7	
Red (640)	660/20	APC	AF647
	780/60	APC-Cy7	APC-eF780